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TA-SE102-AN-KNXX02

096943

Instructions for Use of the SE102-AN pH Sensor

▲ WARNING – Failure to observe this warning may result in serious injury.

Read these instructions for use, observe the Specifications, and follow the Safety Instructions.

1 Safety Instructions

Hazards due to temperature or aggressive media are possible, depending on the location of use. Therefore, the operation and servicing of the sensor shall only be carried out by suitably trained personnel authorized by the operating company. Be sure to wear safety goggles and gloves when cleaning or calibrating a sensor using acid or alkaline solutions.

2 Sensor Description and Intended Use

The SE102-AN pH sensor is equipped with 1 m fixed cable. Connection to the transmitter (e.g. Portavo) is made via a DIN 19262 coax connector and a 4 mm spring connector. The integrated temperature detector provides automatic temperature compensation during calibration and operation of the sensor. The sensor is designed for simultaneous measurement of pH and temperature in industrial processes, particularly in the field of: drinking water, wastewater, water processing. The markings on the sensor or on the packaging label include the following information:

Knick >	Manufacturer of sensor
SE102-AN	Model
0 ... 14	pH range
-5 ... 100 °C	Temperature range
CE	CE conformity mark

The quality certificate and the packaging label provide the serial number and the sensor data of the final inspection. The sensor fits into the sensor quiver of a portable transmitter from the Portavo series. The O-ring (3) seals the quiver.

3 Commissioning

- On unpacking, check the sensor for mechanical damage. Report any damage to your Knick service team.
- Remove the watering cap and briefly rinse the sensor with pure water. After rinsing, the sensor should only be dabbed dry. Do not rub the pH-sensitive glass, since this can lead to electrostatic charging and sluggish response times.
- Connect the sensor cable with the transmitter.

4 Operation

Note: Always open the closure (4) before calibrating or measuring. Refill electrolyte if necessary.

4.1 Calibrating the Sensor

2-point calibration is recommended for the SE102-MS sensor. The sensor is dipped successively into two different buffer solutions with given pH values (e.g., pH 7.00 and pH 4.00). Please refer to the user manual of the pH transmitter for further details.

4.2 Positioning

The sensor can be operated vertically or inclined up to 45° off horizontal.

5 Maintenance and Cleaning

Never keep or store the sensor dry! You should store the sensor with fitted watering cap containing electrolyte (3 mol/l KCl). If the sensor is to be inserted into the quiver of a Portavo series measuring device, the screwed-on watering cap (9) must be replaced by the included flexible watering cap (10). Close the filling hole (4) when storing the sensor. Clean the sensor (glass membrane and junction) with clean water after each production cycle and regularly in case of heavy pollution.

6 Specifications

Model:	SE102-AN
pH:	0 ... 14
Temperature:	-5 ... 100 °C
Length:	110 mm
Junction:	Ceramic
Reference system:	Ag/AgCl with liquid electrolyte

Electrolyte, refillable: 3 mol/l KCl

7 Recommended pH Buffer and Electrolyte Solutions


		Order No.
CaliMat buffer	4.00 (20 °C)	CS-P0400/250
CaliMat buffer	7.00 (20 °C)	CS-P0700/250
Electrolyte	3 mol/l KCl	ZU 0960

For more buffer solutions and electrolytes, see www.knick.de

8 Disposal

Observe the applicable local or national regulations for disposal.


9 Connector plug

Knick >

Manual **SE102-AN**

Read before commissioning. Keep for future use.

- ① Sensor head
- ② Serial number
- ③ O-ring seal for quiver
- ④ **Filling hole with closure**
- ⑤ Nameplate
- ⑥ Glass body
- ⑦ Junctions
- ⑧ pH glass
- ⑨ Watering cap, transport lock
- ⑩ "Watering cap, for use with Portavo quiver"

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Copyright 2020 • Subject to change
 Version: 2 • Published on: March. 03, 2020

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Knick Elektronische Messgeräte GmbH & Co. KG