Portavo® 902 PH

Quickstart guide3

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Instructions succinctes..39

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Guida rapida..... 75

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Version: 1.0

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Web: www.knick.de

Warranty

Defects occurring within 3 years from delivery date shall be remedied free of charge at our plant (carriage and insurance paid by sender).

Sensors and accessories: 1 year

Subject to change

Return of products under warranty

Please contact our Service Team before returning a defective device.

Ship the cleaned device to the address you have been given.

If the device has been in contact with process fluids, it must be decontaminated/ disinfected before shipment. In that case, please attach a corresponding certificate, for the health and safety of our service personnel.

Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

Registered trademarks

The following names are registered trademarks. For practical reasons they are shown without trademark symbol in this manual.

- Calimatic®
- Memosens®
- Paraly®
- Portavo®
- Sensocheck®
- Sensoface®

Check the shipment for transport damage and completeness (see Package Contents).

Caution!

Do not operate the device when one of the following conditions applies:

- the device shows visible damage
- the device fails to perform the intended function
- prolonged storage at temperatures above 70 °C
- severe transport stresses

In this case, a professional routine test must be performed.

This test should be carried out by the manufacturer.

Inserting the Batteries

With four AA batteries, the Portavo has an operating time of over 1000 h.

Open the battery compartment on the rear of the device. Be sure to observe the correct polarity when inserting the batteries (see markings in the battery chamber). Close the battery compartment cover and screw it handtight.

A battery icon in the display indicates the battery power level:

Icon fully filled

Batteries at full capacity

Icon partially filled

Battery capacity is sufficient

Icon empty

Battery capacity not sufficient;
calibration is possible

Icon blinks

Max. 10 operating hours remaining,
measurement is still possible

Caution! It is absolutely necessary to replace the batteries.

Connecting a Sensor

The Portavo 902 PH provides several connections so that many types of sensors can be used for measurement. Note that only **one** sensor may be connected to the meter at a time. The meter automatically recognizes a connected Memosens sensor and switches accordingly. Memosens is signaled in the display.

Separate temperature probe


After power-on, a separate temperature probe is automatically recognized. When you want to replace the temperature probe, you must switch off the meter and then switch it on again.

Caution!

Always make sure that a sensor is connected to the meter before starting measurement.

Explanation: The analog pH input of the Portavo is an electrometer amplifier with an extremely high-impedance. When the sensor is not in contact with the medium or not connected to the meter, electric charges on the input can generate arbitrary, stable pH or mV values which will be shown in the display.


Connections

- a - M8, 4 pins for Memosens sensors
- b - Temperature probe GND
- c - Temperature probe
- d - pH socket (DIN 19 262)


Memosens sensors have a **cable coupling**, which allows convenient replacement of sensors while the cable remains connected to the meter. The connecting cable is connected to socket **a** (M8, 4 pins for Memosens sensors).

Switching On the Meter

When you have connected the sensor, you can switch the meter on by pressing the **on/off** or **meas** key.


When the meter is switched on with the **on/off** key, first a self test is performed and then the calibration data and settings are displayed before the meter switches to measuring mode.


When the meter is switched on with the **meas** key, it immediately switches to measuring mode.

Depending on the connected sensor and the specific measuring task, several steps for configuration and calibration must be performed as described on the following pages.

Icons

Important information about the state of the device:

Keypad

The keys of the membrane keypad have a noticeable pressure point.

They have the following functions:

- on/off** Switches the meter on and displays the device and calibration data (see Start-Up)
- meas** Switches the meter on / Activates measuring mode
- cal** Starts calibration
- set** Activates configuration/ Confirms entries
- clock** Displays time and date, allows setting the clock using **set**
- ▲▼** When this icon is displayed, you can use the arrow keys for navigation.

pH Configuration

Prior to measurement, a configuration should be performed to match the connected sensor and the desired measurement performance.

Furthermore, you can select the suitable calibration method.

The following table gives you an overview. Factory settings are shown in **bold print**.

Measurement

↓ **set**

“Setup” display

Display 1
Display 2
CAL Timer
CAL
CAL POINTS
▲ BUFFER SET (CALIMATIC, FREE CAL) ▼
Auto OFF
Temp Unit
Time Format
Date Format
Default

set

Select using arrow keys, confirm by pressing **set**.

pH x.xx / pH x.xxx / mV / (°C for analog pH only)
OFF / date + time / date / time
OFF / 1 ... 99 days
CALIMATIC /Manual/DATA INPUT/(ISFET-Zero)/FREE CAL
1 / 2 / 3 / 1-2-3 (for CALIMATIC, Manual, FREE CAL)
-01- Mettler Toledo 2.00/4.01/7.00/9.21
-02- Knick CaliMat 2.00/4.00/7.00/9.00/12.00
-03- Ciba (94) 2.06/4.00/7.00/10.00
-04- NIST technical 1.68/4.00/7.00/10.01/12.46
-05- NIST standard 1.679/4.006/6.865/9.180
-06- HACH 4.01/7.00/10.01/12.00
-07- WTW techn. buffers 2.00/4.01/7.00/10.00
-08- Hamilton 2.00/4.01/7.00/10.01/12.00
-09- Reagecon 2.00/4.00/7.00/9.00/12.00
-10- DIN 19267 1.09/4.65/6.79/9.23/12.75
OFF / 0.1h / 1h / 6h / 12h
°C / °F
24h / 12h
dd.mm.yy / mm.dd.yy
NO / YES (reset to factory settings)


- ▲ This icon prompts you to select a menu item using the arrow keys –
- ▼ the selection is confirmed by pressing **set**.

CALIMATIC Calibration

(Calibration with automatic buffer recognition)

The calibration method is selected in the configuration menu. Calibration is required to adjust the sensor to the meter. It is indispensable for achieving comparable and reproducible measurement results.


Note: To abort calibration, you can press **meas** at any time. This will be confirmed by the display message "CAL ABORTED". Exception: When you have selected "CAL POINTS 1-2-3" and the first calibration step has been completed, the calibration process cannot be stopped any more.

DATA INPUT Calibration

(Calibration by entering known sensor values)

The calibration method is selected in the configuration menu.

The calibration data will be displayed successively:

Date and time

ZERO POINT

SLOPE

Then the meter switches to measuring mode.

Note: To abort calibration, you can press **meas** at any time.

MANUAL Calibration

(Manual calibration)

The calibration method is selected in the configuration menu.

Measurement

↓ **cal**

CAL
MANUAL

The number of calibration points has been selected in the configuration menu.

↓

CAL 1/2/3
PRESS CAL

↓ **cal**

pH display blinks
PRESS CAL

Use ▲▼ to set the temperature-corrected pH value taken from the buffer table.

↓ **cal**

mV display blinks

Depending on the number of calibration points, the procedure described above for CAL 1/2/3 is repeated.

↓

mV value blinks until calibration is completed, then successive display of:

CAL DATA

1/2/3 CAL POINTS

ZERO POINT

SLOPE

Then the meter switches to measuring mode.

Note: To abort calibration, you can press **meas** at any time. This will be confirmed by the display message "CAL ABORTED". Exception: When you have selected "CAL POINTS 1-2-3" and the first calibration step has been completed, the calibration process cannot be stopped any more.

FREE CAL Calibration

(Free selection of calibration method)

FREE CAL calibration is selected in the configuration menu.

Measurement

↓ cal


CAL
CALIMATIC blinks

↓ cal


Use ▲▼ to select the required calibration method (CALIMATIC, DATA INPUT or MANUAL).

Perform the selected calibration (see CALIMATIC, DATA INPUT or MANUAL calibration).

Error messages are indicated as “ERROR ...” on the display. Information on the sensor condition is indicated by the “Sensoface” icon (friendly, neutral, sad) possibly accompanied by an info message (“INFO ...”).

Example of an error message:
ERROR 8 (identical calibration media)

Example of a “Sensoface” message:
INFO 1 (cal timer expired)

Sensoface (the “smiley” icon) provides information on the sensor condition (maintenance request). Measurement can still be performed. After a calibration, the corresponding Sensoface icon (friendly, neutral, sad) is shown together with the calibration data. Otherwise, Sensoface is only visible in measuring mode.

The most important error messages and “Sensoface” info messages are shown on the inside of the protective cover.

A complete list of messages and their meanings is provided in the following tables.

“Sensoface” Messages

The “Sensoface” icon provides information on the sensor condition:

Sensoface Meaning

Sensor is okay

Calibrate the sensor soon

Calibrate or replace the sensor

The “neutral” and “sad” Sensoface icons are accompanied by an “INFO ...” message to give a hint to the cause of deterioration.

Sensoface	Message	Cause
	INFO 1	Calibration timer
	INFO 3	Sensocheck
	INFO 5	Zero / Slope
 	INFO 6	Response time
	INFO 7	ISFET: Operating point (asymmetry potential)
	INFO 8	ISFET: Leakage current
	INFO 9	ORP offset

Error Messages

The following error messages can be shown in the display.

 blinks	Battery empty	Replace batteries
ERROR 1	pH value out of range	Check whether the measurement conditions correspond to the adjusted measuring range.
ERROR 2	ORP value out of range	
ERROR 3	Temperature value out of range	
ERROR 4	Sensor zero point too high/low	Thoroughly rinse the sensor and recalibrate. If this does not help, replace the sensor.
ERROR 5	Sensor slope too high/low	
ERROR 8	Calibration error: Identical buffers	Use a buffer solution with a different nominal value before starting the next calibration step.
ERROR 9	Calibration error: Buffer unknown	Make sure that you use the same buffer set as configured.
ERROR 10	Cal media interchanged	Repeat calibration.
ERROR 11	Measured value unstable Drift too high	Leave the sensor in the liquid until the temperature is stable. If this does not help, replace the sensor.
ERROR 14	Time and date invalid	Set time and date
ERROR 18	Configuration invalid	Restart, reset to factory settings (Setup: DEFAULT YES), configure and calibrate. If this does not help, send in the device for repair.
ERROR 19	Factory settings error	Device defective, send it in.
ERROR 21	Sensor error (Memosens)	Connect operational Memosens sensor.

pH/mV input	pH socket, DIN 19 262 (13/4 mm)	
pH range	-2 ... 16	
Decimal places ^{*)}	2 or 3	
	Input resistance	1 x 10 ¹² Ω (0 ... 35 °C)
	Input current	1 x 10 ⁻¹² A (at RT, doubles every 10 K)
Measuring cycle	Approx. 1 s	
Measurement error ^{1,2,3)}	< 0.01 pH, TC < 0.001 pH/K	
mV range	-1300 ... +1300 mV	
Measuring cycle	Approx. 1 s	
Measurement error ^{1,2,3)}	< 0.1 % meas. val. + 0.3 mV, TC < 0.03 mV/K	
Temperature input	2 x 4 mm dia. for integrated or separate temperature detector	
Measuring ranges	NTC30 temp detector	-20 ... +120 °C
	Pt1000 temp detector	-40 ... +250 °C
Measuring cycle	Approx. 1 s	
Measurement error ^{1,2,3)}	< 0.2 K (Tamb = 23 °C); TC < 25 ppm/K	
Memosens pH input	M8 socket, 4 pins, for Memosens lab cable	
Display ranges ⁴⁾	pH	-2.00 ... +16.00
	mV	-2000 ... +2000 mV
	Temperature	-50 ... +250 °C
Memosens pH input ISFET	M8 socket, 4 pins, for Memosens lab cable	
Display ranges ⁴⁾	pH	-2.00 ... +16.00
	mV	-2000 ... +2000 mV
	Temperature	-50 ... +250 °C
Memosens ORP input	M8 socket, 4 pins, for Memosens lab cable	
Display ranges ⁴⁾	mV	-2000 ... +2000 mV
	Temperature	-50 ... +250 °C
Sensor standardization ^{*)}	ORP calibration (zero adjustment)	
Permissible calibration range	ΔmV (offset)	-700 ... +700 mV

*) User-defined

1) According to EN 60746-1,
at nominal operating conditions
2) ± 1 count

3) Plus sensor error
4) Ranges depending on Memosens sensor

Sensor standardization *)	pH calibration	
Operating modes *)	CALIMATIC	Calibration with automatic buffer recognition
	MANUAL	Manual calibration with entry of individual buffer values
	DATA INPUT	Data entry of zero and slope
Calimatic buffer sets *)	-01- Mettler-Toledo	2.00/4.01/7.00/9.21
	-02- Knick CaliMat	2.00/4.00/7.00/9.00/12.00
	-03- Ciba (94)	2.06/4.00/7.00/10.00
	-04- NIST technical	1.68/4.00/7.00/10.01/12.46
	-05- NIST standard	1.679/4.006/6.865/9.180
	-06- HACH	4.01/7.00/10.01 /12.00
	-07- WTW techn. buffers	2.00/4.01/7.00/10.00
	-08- Hamilton	2.00/4.01/7.00/10.01/12.00
	-09- Reagecon	2.00/4.00/7.00/9.00/12.00
	-10- DIN 19267	1.09/4.65/6.79/9.23/12.75
Permissible calibration range	Zero point	pH 6 ... 8
	With ISFET:	-750 ... +750 mV
	Operating point (asymmetry)	
	Slope	approx. 74 ... 104 % (possibly restricting notes from Sensoface)
Calibration timer *)	Interval 1 ... 99 days, can be switched off	
Sensoface	Provides information on the sensor condition	
Evaluation of	zero/slope, response, calibration interval	

*) User-defined

Connections	1 x pH socket, DIN 19 262 2 x 4-mm socket for separate temperature detector 1 x M8 socket, 4 pins, for Memosens lab cable
Display	LCD STN 7-segment display with 3 lines and icons
Sensoface	Status indication (friendly, neutral, sad)
Status indicators	Battery power level
Notices	Hourglass
Keypad	[on/off], [cal], [meas], [set], [▲], [▼], [clock]
Diagnostics functions	
Sensor data (Memosens only)	Manufacturer, sensor type, serial number, operating time
Calibration data	Calibration date, zero, slope
Device self-test	Automatic memory test (FLASH, EEPROM, RAM)
Device data	Device type, software version, hardware version
Data retention	Parameters, calibration data > 10 years
EMC	EN 61326-1 (General Requirements)
Emitted interference	Class B (residential area)
Immunity to interference	Industry EN 61326-2-3 (Particular Requirements for Transmitters)
RoHS conformity	According to directive 2011/65/EC
Power supply	
Portavo 902	4 x AA alkaline batteries
Operating time	Approx. 1000 h (alkaline)
Nominal operating conditions	
Ambient temperature	-10 ... +55 °C
Transport/ Storage temperature	-25 ... +70 °C
Relative humidity	0 ... 95 %, short-term condensing allowed
Housing	
Material	PA12 GF30 (silver gray RAL 7001) + TPE (black)
Protection	IP 66/67 with pressure compensation
Dimensions	Approx. (132 x 156 x 30) mm
Weight	Approx. 500 g

A

AA batteries 4
Analog pH input 5
Arrow keys 7
Automatic calibration (Calimatic) 9

B

Battery capacity 4
Battery charge indicator 4
Battery compartment 4
Battery icon 4
Battery replacement 4
Buttons 7

C

Calibration (pH), Calimatic 9
Calibration (pH), data input 10
Calibration (pH), manual 11
Calibration (pH), method selected in measuring mode 12
Calimatic automatic calibration 9
cal key 7
Charge level of batteries 4
clock key 7
Configuration (pH) 8
Connecting a sensor 5
Connecting cable for Memosens 5
Connections 5
Control elements 7

D

Data input (pH calibration) 10
Data of the meter 16
Device configuration 8
Device messages 13
Display icons 6
Disposal 3

E

ERROR (error codes) 15
Error messages 13
Error messages, overview 15

F

FREE CAL 12

I

Icons in display 6
INFO messages 14
Inserting the batteries 4

K

Keypad 7

M

Manual calibration 11
meas key 7
Memosens connecting cable 5
Menu structure of configuration 8
Messages 13

O

on/off key 7
Overview of configuration 8
Overview of error messages 15

P

Parameter settings (configuration) 8
pH configuration 8
pH input, analog 5
pH socket, DIN 19 262 5
Ports 5
Power-on 6

R

Registered trademarks 3
Replacing the batteries 4
Return of products under warranty 3

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© 2017 Änderungen vorbehalten