Portavo® 907 MULTI

Quickstart guide3

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Instructions succinctes 109

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Other languages: www.knick.de

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EN Copyright 2019 • Subject to change

Version: 4

This document was last updated on February 7, 2019

The latest documents are available for download on our website under the corresponding product description.

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TA-209.7MU-KNX04

Software version: 1.x

Start-Up	6
Inserting the Batteries	6
Connecting a Sensor	7
Switching On the Meter	8
Icons	8
Display and Keypad	9
Overview of the Portavo 907	9
Information	10
Calibration Record	10
Sensor Information (Digital Sensors only)	10
Sensor Network Diagram (pH and Oxy only)	11
Sensor Monitor	12
Messages	12
MemoLog (Memosens only)	12
Device Info	13
Device Test	13
pH Configuration	14
ORP Configuration	16
Conductivity Configuration	18
Oxygen Configuration	21
pH Calibration	23
Calimatic Calibration	23
Manual Calibration	24
Data Entry Calibration	24
Calibrating a pH/ORP Combo Sensor	24
ORP Calibration	25
ISFET Calibration	26
Conductivity Calibration	27
Auto Calibration	27
“Entry of Solution” Calibration	28
Cell Constant / Cell Factor Calibration	28
Installation Factor Calibration	28
Zero Calibration	28
Oxygen Calibration	29
Calibration in Air	29
Zero Calibration	30
Data Entry Calibration	30

Measuring	31
Switching the Measured Value Display	31
Adjusting the Temperature.....	31
Data Logger	32
Operating Modes of the Data Logger (Logger Type)	33
Configuring the Data Logger.....	35
Increasing the Battery Life	35
Starting/Stopping the Data Logger.....	37
Viewing the Logger Data.....	37
Deleting the Logger Data.....	37
Error and Status Messages	39
“Sensoface” Messages.....	40
Info and Help Texts	40
Error Messages	41
Option 001 SOP	42
Option 002 Temp.Cal	42
How to Enter the Rescue PIN Code	46
Specifications	48
Index	54

Return of Products Under Warranty

Please contact our Service Team before returning a defective device.

Ship the cleaned device to the address you have been given.

If the device has been in contact with process fluids, it must be decontaminated/ disinfected before shipment. In that case, please attach a corresponding certificate, for the health and safety of our service personnel.

Disposal

Please observe the applicable local or national regulations concerning the disposal of "waste electrical and electronic equipment".

Registered Trademarks

The following names are registered trademarks. For practical reasons they are shown without trademark symbol in this manual.

- CaliMat®
- Calimatic®
- Memosens®
- Paraly®
- Portavo®
- Sensocheck®
- Sensoface®

Check the shipment for transport damage and completeness (see Package Contents).

NOTICE!

Do not operate the device when one of the following conditions applies:

- the device shows visible damage
- the device fails to perform the intended function
- prolonged storage at temperatures above +70 °C / +158 °F
- severe transport stresses

In this case, a professional routine test must be performed.

This test should be carried out at our factory.

Inserting the Batteries

With four AA batteries, the Portavo has an operating time of up to 500 h when operated in logger mode (see page 35).

Open the battery compartment on the rear of the device. Be sure to observe the correct polarity when inserting the batteries (see markings in the battery chamber). Close the battery compartment cover and screw it handtight.

A special lithium-ion battery (ZU 0925) suited to the battery compartment is available for the Portavo 907. Only this battery type can be charged directly from the USB port.

A battery icon in the display indicates the battery power level:

	Icon fully filled	Batteries at full capacity
	Icon partially filled	Battery capacity is sufficient
	Icon empty	Battery capacity not sufficient; calibration is possible, no logging
	Icon blinks	Only a few operating hours remaining, measurement is still possible NOTICE! It is absolutely necessary to replace the batteries.

Connecting a Sensor


The Portavo 907 MULTI provides a special socket for connecting conventional sensors for the respective process parameter. Alternatively, you can connect a Memosens sensor for pH, conductivity or oxygen measurement. The Model SE 340 optical oxygen sensor can also be connected. The meter automatically recognizes a connected Memosens sensor and accordingly selects the corresponding process variable. Memosens is signaled in the display.

Note that only **one** sensor may be connected to the meter at a time.

Separate temperature probe

After power-on, a separate temperature probe is automatically recognized.

When you want to replace the temperature probe, you must switch off the meter and then switch it on again.


Connections

- a - Micro USB port
- b - M8, 4 pins, for Memosens lab cable
- c - Temperature probe GND
- d - Temperature probe
- e - Depending on the device variant:
 - Portavo 907 MULTI PH: pH socket acc. to DIN 19 262 for analog sensors
 - Portavo 907 MULTI COND: DIN socket, 8 pins, for analog sensors
 - Portavo 907 MULTI OXY: M12, 8 pins, for Memosens sensors or SE 340 sensor (optical oxygen)

Memosens sensors have a **cable coupling** which allows convenient replacement of sensors while the cable remains connected to the meter. The connecting cable is connected to socket **b** (Memosens lab cable) or **e** (flexible connecting cable – for Portavo 907 MULTI OXY only!).

Switching On the Meter

You can use **meas** or **on/off** to switch the meter on:

Analog sensors:

- When you press **meas** or **on/off**, the meter immediately switches to measuring mode.


Memosens sensors:

- When you press **meas** or **on/off**, the meter displays selected sensor data before it switches to measuring mode.


Icons

Important information about the state of the device:

Display and Keypad

Display and keypad correspond directly via softkeys.

Menu icons

Information

Calibration

Data logger

Configuration

- Softkeys** Function is shown in the display above the key
- Arrow keys** Selecting / Adjusting entries
- enter** Confirming an adjustment
- on/off** Switching on / off
- meas** Switching on / Immediate access of meas. mode /
Toggling the display / Displaying time and date

pH

ORP

Oxy

Cond

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Information" and confirm by pressing **enter**.
- 3) Select the desired submenu and confirm by pressing **enter**.

The different submenus are described below.

Calibration Record

Shows the data of the last calibration performed on the currently connected sensor.

Sensor Information (Digital Sensors only)

Shows the data of the currently connected sensor. When MemoLog has been activated (in the Configuration menu), you can save the sensor data in the device by pressing the **Save** softkey. The following table shows the sensor information depending on the respective sensor:

	pH/ pH/ORP**	Cond	Oxy	ISFET	ORP	Optical Oxy
Manufacturer	x	x	x	x	x	x
Ref. No.	x	x	x	x	x	x
Sensor serial no.	x	x	x	x	x	x
Cap serial no.						x
TAG	x	x	x	x	x	
SW version	x	x	x	x	x	x
HW version	x	x	x	x	x	
Calibration*	x	x	x	x	x	x
Zero point	x		x			x
Slope	x		x	x		x
ORP calibration***	x					
Correction					x	
Nom. cell constant		x				
Temp. offset	x	x	x		x	
Sensor operating time	x	x	x	x	x	x
Cap operating time						x
Wear	x		x	x		
SIP	x	x	x	x	x	
CIP	x**	x				
Autoclaving	x**					
Cell constant		x				
Operating point				x		

* latest calibration ** for pH/ORP combo sensor only

pH

Oxy

Sensor Network Diagram (pH and Oxy only)


Provide single-glance information on the following parameters of the connected sensor:

- Slope
- Zero point (operating point for Memosens ISFET)
- Sensocheck (pH) or leakage current (ISFET and Oxy)
- Response time
- Calibration timer
- Wear (Memosens)

Parameters which cannot be checked are shown as inactive (gray) and are set to 100% (e.g. Sensocheck for analog sensors).

The parameter values should lie between the outer (100%) and inner (50%) hexagon. When a value enters the inner hexagon (<50%), the corresponding caption text flashes red (see example).

Example: Network diagram of a digital pH sensor (Memosens)

pH

ORP

Oxy

Cond

Sensor Monitor

Shows the raw values available from the connected sensor:

pH, analog	mV, temperature, temperature detector, temp. resistance
pH, digital, glass	mV, temperature, glass impedance
pH, digital, ISFET	mV, leakage current, temperature
pH, ORP	mV, temperature
Cond, analog	Resistance, conductance, temperature, temperature detector, temperature resistance
Cond, digital	Resistance, conductance, temperature
Oxy, digital	Sensor current, leakage current, polarization voltage, partial pressure, air pressure, temperature
Oxy, digital, optical	Partial pressure, temperature

Messages

Shows all active error and status messages as well as supplementary help texts.

MemoLog (Memosens only)

Displays the individual calibration records stored in the device. You have the possibility to delete individual entries or all entries. The following parameters are displayed:

- Sensor type
- Serial no.
- TAG
- Calibration date
- Zero point
- Slope
- Cell constant (Cond sensor)
- Operating point (ISFET sensor)

Background: The device provides a calibration data logger, which must be activated in the configuration menu. With "MemoLog" activated, up to 100 calibration records can be directly saved to the device. After every calibration, the complete Memosens data will be recorded. Convenient management of the calibration data is possible using the MemoSuite or Paraly SW 112 software.

MemoLog is not suitable for SE 340 (optical oxygen sensor).

pH

ORP

Oxy

Cond

Device Info

Shows the following device information:

- Model name
- Serial number
- Software version
- Hardware version
- Air pressure
- Battery

Device Test

A device self-test is automatically run in the background at regular intervals.

It checks the memory modules listed below. A green checkmark shows that the test was successful.

- FLASH program memory
- FLASH data memory
- FLASH parameter memory
- RAM (working memory)

Display test

- 1) Select "Display test" and press **enter**.
- 2) The display lights up red, green, blue and then white.
- 3) Press any key to stop the test.

Keypad test

- 1) Select "Keypad test" and press **enter**.
- 2) Press all nine keys one after the other.
A green checkmark shows that a key functions properly.
- 3) Press any key to stop the test.

pH

pH Configuration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Make the desired adjustments.

The following table gives you an overview.

Factory settings are shown in **bold print**.

"pH Configuration" menu selection – part 1

Language	Deutsch English Español Italiano Français Português Русский
Auto-off	Off 5 min 10 min 30 min 60 min
Temperature	°C °F
Right softkey	Logger Start/Stop Hold value
+ pH sensor*	
Display format	0.00 pH 0.000 pH
Wear	On Off
+ Calibration*	
Cal timer	Off On
Interval	On: 00 ... 99 days
Cal mode	Calimatic Manual Data entry
Cal points	Auto 1-point 2-point 3-point
Buffer set	Mettler-Toledo 2.00/4.01/7.00/9.21 Knick CaliMat 2.00/4.00/7.00/9.00/12.00 Ciba 2.06/4.00/7.00/10.00 NIST technical 1.68/4.00/7.00/10.01/12.46 NIST standard 1.679/4.006/6.865/9.180 Hach 4.01/7.00/10.01/12.00 WTW 2.00/4.01/7.00/10.00 Hamilton 2.00/4.01/7.00/10.01/12.00 Reagecon 2.00/4.00/7.00/9.00/12.00 DIN 19267 1.09/4.65/6.79/9.23/12.75 Metrohm 4.00/7.00/9.00 User buffer 1**
MemoLog	Off On
TAG	Off On

* "+" indicates that submenus can be opened by pressing **enter**.

** Parameter can be configured using Paraly SW 112 software.

“pH Configuration” menu selection – part 2

	+ Time/Date*	
	Time format	24 h 12 h
	Date format	dd.mm.yyyy yyyy-mm-dd dd/mm/yyyy mm/dd/yyyy
	Time	hh:mm:ss
	Date	Date format as configured
	+ Display*	
	Appearance	Modern Retro
	Lighting	Permanent 60 min 30 min 10 min 5 min 1 min 30 sec
	Brightness	Bright Standard Dim
	+ Data logger*	
	Meas.point	--
	Note	--
	Recording	Non-circular Circular
	Logger type	Shot
		Interval 00.00.01...12:59:59 00:02:00
	Difference 1st difference On Off	
	Delta pH pH 0.0...16.0 pH 1.0	
	Delta mV 0 ... 2000 mV 1 mV	
	2nd diff. On Off	
	Delta °C 0...99.9 °C 1.0 °C	
	Delta °F 0...450 °F 1.0 °F	
	Intv+Diff Interval See logger type: Interval	
	Difference See logger type: Difference	
	Limit value Interval Basis/Event 00.00.01...12:59:59 00:01:00/00:00:01	
	Limit values Min/Max, corresponding to permissible range (see Specifications)	
+ Options	001 SOP Add-on function, enabled 002 Temp.cal via TAN	
Factory setting	Yes No Note: Reset to factory settings will also erase all logger data!	

* “+” indicates that submenus can be opened by pressing **enter**.

ORP

ORP Configuration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Make the desired adjustments.

The following table gives you an overview.

Factory settings are shown in **bold print**.

"ORP Configuration" menu selection – part 1

Language	Deutsch English Español Italiano Français Português Русский
Auto-off	Off 5 min 10 min 30 min 60 min
Temperature	°C °F
Right softkey	Logger Start/Stop Hold value
+ Calibration*	
MemoLog	Off On
TAG	Off On
+ Time/Date*	
Time format	24 h 12 h
Date format	dd.mm.yyyy yyyy-mm-dd dd/mm/yyyy mm/dd/yyyy
Time	hh:mm:ss
Date	Date format as configured
+ Display*	
Appearance	Modern Retro
Lighting	Permanent 60 min 30 min 10 min 5 min 1 min 30 sec
Brightness	Bright Standard Dim

* "+" indicates that submenus can be opened by pressing **enter**.

Cond

Conductivity Configuration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Make the desired adjustments.

The following table gives you an overview.


Factory settings are shown in **bold print**.

"Conductivity Configuration" menu selection – part 1

Language	Deutsch English Español Italiano Français Português Русский
Auto-off	Off 5 min 10 min 30 min 60 min
Temperature	°C °F
Right softkey	Logger Start/Stop Hold value
+ Cond sensor*	
Conductivity	S/cm S/m
Range selection	Auto 0.000 µS/cm 00.00 µS/cm 000.0 µS/cm 0000 µS/cm 00.00 mS/cm 000.0 mS/cm 0000 mS/cm
Compensation	Off MΩ cm TC SAL TDS % by wt
TC comp.	TC: Linear NLF NaCl HCl NH3 NaOH
TC of solution	TC: 0 ... 20.0 %/K 2.1 %/K
Ref. temp	TC: 0 ... 100.0 °C 25 °C 32 ... 212 °F 77 °F
TDS factor	TDS: 0 ... 9.99 1.00
Solution	% by wt: NaCl HCl NaOH H2SO4 HNO3

* "+" indicates that submenus can be opened by pressing **enter**.

“Conductivity Configuration” menu selection – part 2

* "+" indicates that submenus can be opened by pressing **enter**.

Cond

“Conductivity Configuration” menu selection – part 3

<div style="display: flex; align-items: center;"> <div style="margin-right: 5px;">↑</div> <div style="border: 1px solid black; padding: 2px;">+ Data logger*</div> </div> <div style="border: 1px solid black; padding: 2px; margin-top: 2px;">Meas.point</div> <div style="border: 1px solid black; padding: 2px; margin-top: 2px;">Note</div> <div style="border: 1px solid black; padding: 2px; margin-top: 2px;">Recording</div> <div style="border: 1px solid black; padding: 2px; margin-top: 2px;">Logger type</div> <div style="border: 1px solid black; padding: 2px; margin-top: 2px; height: 100px;"></div> <div style="display: flex; align-items: center; margin-top: 2px;"> <div style="margin-right: 5px;">↑</div> <div style="margin-right: 5px;">▼</div> </div> <div style="border: 1px solid black; padding: 2px; margin-top: 2px;">+ Options</div> <div style="border: 1px solid black; padding: 2px; margin-top: 2px;">Factory setting</div> <div style="margin-top: 5px;">↓</div>	<div style="display: flex; flex-direction: column; align-items: center;"> <div style="margin-bottom: 5px;">▲</div> <div style="margin-bottom: 5px;">▼</div> <div style="margin-bottom: 5px;">↻</div> <div style="margin-bottom: 5px;">↔</div> </div> <p>enter</p>	<table border="1" style="border-collapse: collapse; width: 100%;"> <tr><td colspan="3">--</td></tr> <tr><td colspan="3">--</td></tr> <tr><td colspan="3">Non-circular Circular</td></tr> <tr><td colspan="3">Shot</td></tr> <tr><td colspan="3">Interval 1...12:59:59 00:02:00</td></tr> <tr> <td rowspan="6" style="width: 15%;">Difference</td> <td style="width: 30%;">1st difference</td> <td>On Off</td> </tr> <tr> <td>Delta cond</td> <td>0 ... 2000 mS/cm 1.0 μS/cm</td> </tr> <tr> <td>Delta conc</td> <td>0 ... 9.99 % 1.0 %</td> </tr> <tr> <td>Delta MΩcm</td> <td>0 ... 9.999 MΩcm 1 MΩcm</td> </tr> <tr> <td>Delta salinity</td> <td>0 ... 45.0 g/kg 1.0 g/kg</td> </tr> <tr> <td>Delta TDS</td> <td>0 ... 2000.0 mg/l 1 mg/l</td> </tr> <tr> <td></td> <td>2nd diff.</td> <td>On Off</td> </tr> <tr> <td></td> <td>Delta °C</td> <td>0...99.9 °C 1.0 °C</td> </tr> <tr> <td></td> <td>Delta °F</td> <td>0...450 °F 1.0 °F</td> </tr> <tr> <td rowspan="2">Intv+Diff</td> <td>Interval</td> <td>See logger type: Interval</td> </tr> <tr> <td>Difference</td> <td>See logger type: Difference</td> </tr> <tr> <td rowspan="2">Limit value</td> <td>Interval</td> <td>Basis/Event 00.00.01...12:59:59 00:01:00/00:00:01</td> </tr> <tr> <td>Limit values</td> <td>Min/Max, corresponding to permissible range (see Specifications)</td> </tr> <tr> <td colspan="2">001 SOP</td> <td>Add-on function, enabled via TAN</td> </tr> <tr> <td colspan="2">002 Temp.cal</td> <td></td> </tr> <tr> <td colspan="3">Yes No</td> </tr> <tr> <td colspan="3">Note: Reset to factory settings will also erase all logger data!</td> </tr> </table>	--			--			Non-circular Circular			Shot			Interval 1...12:59:59 00:02:00			Difference	1st difference	On Off	Delta cond	0 ... 2000 mS/cm 1.0 μS/cm	Delta conc	0 ... 9.99 % 1.0 %	Delta MΩcm	0 ... 9.999 MΩcm 1 MΩcm	Delta salinity	0 ... 45.0 g/kg 1.0 g/kg	Delta TDS	0 ... 2000.0 mg/l 1 mg/l		2nd diff.	On Off		Delta °C	0...99.9 °C 1.0 °C		Delta °F	0...450 °F 1.0 °F	Intv+Diff	Interval	See logger type: Interval	Difference	See logger type: Difference	Limit value	Interval	Basis/Event 00.00.01...12:59:59 00:01:00/00:00:01	Limit values	Min/Max, corresponding to permissible range (see Specifications)	001 SOP		Add-on function, enabled via TAN	002 Temp.cal			Yes No			Note: Reset to factory settings will also erase all logger data!		
--																																																													
--																																																													
Non-circular Circular																																																													
Shot																																																													
Interval 1...12:59:59 00:02:00																																																													
Difference	1st difference	On Off																																																											
	Delta cond	0 ... 2000 mS/cm 1.0 μS/cm																																																											
	Delta conc	0 ... 9.99 % 1.0 %																																																											
	Delta MΩcm	0 ... 9.999 MΩcm 1 MΩcm																																																											
	Delta salinity	0 ... 45.0 g/kg 1.0 g/kg																																																											
	Delta TDS	0 ... 2000.0 mg/l 1 mg/l																																																											
	2nd diff.	On Off																																																											
	Delta °C	0...99.9 °C 1.0 °C																																																											
	Delta °F	0...450 °F 1.0 °F																																																											
Intv+Diff	Interval	See logger type: Interval																																																											
	Difference	See logger type: Difference																																																											
Limit value	Interval	Basis/Event 00.00.01...12:59:59 00:01:00/00:00:01																																																											
	Limit values	Min/Max, corresponding to permissible range (see Specifications)																																																											
001 SOP		Add-on function, enabled via TAN																																																											
002 Temp.cal																																																													
Yes No																																																													
Note: Reset to factory settings will also erase all logger data!																																																													

* “+” indicates that submenus can be opened by pressing **enter**.

Oxygen Configuration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Make the desired adjustments.

The following table gives you an overview.

Factory settings are shown in **bold print**.

"Oxygen Configuration" menu selection – part 1

	Language	Deutsch English Español Italiano Français Português Русский
	Auto-off	Off 5 min 10 min 30 min 60 min
	Temperature	°C °F
	Right softkey	Logger Start/Stop Hold value
	+ Oxy sensor*	
	Display	Saturation Concentration Partial pressure
	Salinity	0 ... 45.0 g/kg
	Pressure correc- tion**	Air pressure Manual
	Pressure	Manual: 0 ... 9999 mbar 1013 mbar
	Wear	On Off
	+ Calibration*	
	Cal mode	In air Zero point Data entry
	Calibration timer	Off On
	Interval	On: 0 ... 99 days
	MemoLog	Off On
TAG	Off On	
+ Time/Date*		
Time format	24 h 12 h	
Date format	dd.mm.yyyy yyyy-mm-dd dd/mm/yyyy mm/dd/yyyy	
Time	hh:mm:ss	
Date	Date format as configured	

* "+" indicates that submenus can be opened by pressing **enter**.

** The device provides an internal barometer.

Oxy

"Oxygen Configuration" menu selection – part 2

	+ Display*
	Appearance
	Lighting
	Brightness
	+ Data logger*
	Meas.point
	Note
	Right softkey
	Recording
	Logger type
+ Options	
Factory setting	

Modern Retro
Permanent 60 min 30 min 10 min 5 min 1 min 30 sec
Bright Standard Dim

--		
--		
Logger Start/Stop Hold value		
Non-circular Circular		
Shot		
Interval	00.00.01...12:59:59 00:02:00	
Difference	1st difference	On Off
	Delta saturation	0 ... 200 %Air 1% Air
	Delta conc	0 ... 20 mg/l 1 mg/l
	Delta mbar	0 ... 999.99 mbar 1 mbar
	2nd diff.	On Off
	Delta °C	0...99.9 °C 1.0 °C
	Delta °F	0...450 °F 1.0 °F
Intv+Diff	Interval	See logger type: Interval
	Difference	See logger type: Difference
Limit value	Interval	Basis/Event 00.00.01...12:59:59 00:01:00/00:00:01
	Limit values	Min/Max, corresponding to permissible range (see Specifications)
001 SOP	Add-on function, enabled via TAN	
002 Temp.cal	Add-on function, enabled via TAN	
Yes No		
Note: Reset to factory settings will also erase all logger data!		

* "+" indicates that submenus can be opened by pressing **enter**.

pH Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) The "TAG" menu allows the sensor TAG to be edited.
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 5) Perform the selected calibration as described on the following pages.
Follow the instructions on the display.

Note: Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

Calimatic Calibration

(Automatic calibration with specification of the buffer solution used)

- 1) Select the number of calibration points and the buffer set as shown in the table below and press the **Start** softkey.

Cal points	Auto	1-point	2-point	3-point
Buffer set	Mettler-Toledo		2.00/4.01/7.00/9.21	
	Knick CaliMat		2.00/4.00/7.00/9.00/12.00	
	Ciba		2.06/4.00/7.00/10.00	
	NIST technical		1.68/4.00/7.00/10.01/12.46	
	NIST standard		1.679/4.006/6.865/9.180	
	Hach		4.01/7.00/10.01/12.00	
	WTW		2.00/4.01/7.00/10.00	
	Hamilton		2.00/4.01/7.00/10.01/12.00	
	Reagecon		2.00/4.00/7.00/9.00/12.00	
	DIN 19267		1.09/4.65/6.79/9.23/12.75	
	Metrohm		4.00/7.00/9.00	
	User buffer 1		configurable via Paraly SW 112 software	

- 2) Immerse the sensor in the **1st/2nd/3rd** buffer solution and press **Continue** (repeat this step for each calibration point).
- 3) Finally, the calibration data will be displayed.
You can **Apply** or **Discard** these values.

Note: To abort calibration, you can press **meas** at any time.

pH

Manual Calibration

(Calibration with manual specification of the number of calibration points and the buffer solution)

- 1) Select the number of calibration points and press the **Start** softkey.
- 2) Adjust the temperature-corrected value (see buffer table) for the **1st/2nd/3rd** buffer solution and press **Continue** (repeat this step for each calibration point).

Note: When using sensors without temperature detector, you should adjust the temperature manually before starting calibration (see page 31).

- 3) Finally, the calibration data will be displayed.

You can **Apply** or **Discard** these values.

Data Entry Calibration

(Calibration by entering known sensor values)

- 1) Press the **Start** softkey.
- 2) Enter the known sensor values for zero and slope.
- 3) Finally, you can **Apply** these values or **Cancel** the calibration.

pH

ORP

Calibrating a pH/ORP Combo Sensor

The pH/ORP combo sensor can be calibrated as a pH sensor and/or as an ORP sensor.

pH Calibration

Follow the instructions given for pH calibration, page 23.

ORP Calibration

Follow the instructions given for ORP calibration, page 25.

Note: To abort calibration, you can press **meas** at any time.

ORP Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) The "TAG" menu allows the sensor TAG to be edited.
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 5) Enter the temperature-corrected setpoint of the calibration solution.
- 6) Immerse the sensor in the calibration solution and wait until the measured value is stable.
- 7) **Apply** or **Discard** the ORP setpoint.

Note: Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

Note: To abort calibration, you can press **meas** at any time.

ISFET

ISFET Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) The "TAG" menu allows the sensor TAG to be edited.
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 5) Perform the selected calibration as described on the following pages.
Follow the instructions on the display.

Note: Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

Calibrating the ISFET Zero (Operating Point)

- 1) Select the "ISFET zero" calibration mode for setting the operating point for the first sensor calibration.

Calibration mode	Calimatic
	Manual
	Data entry
	ISFET zero (operating point)

- 2) Press the **Start** softkey.
- 3) Adjust the buffer value if required: default pH 7.00
- 4) Press the **Start** softkey.
- 5) Finally, you can **Apply** or **Discard** the calibration value for the operating point.
When you apply the calibration value, the operating point will be stored in the device, but not in the sensor!
Keep the sensor connected to the Portavo while performing the next calibration step. The operating point will be taken into account for the following calibration.

Calimatic/Manual/Data Entry Calibration

Follow the instructions given for pH calibration, page 23.

If you disconnect the sensor before performing the calibration (e.g., Calimatic), you must set the operating point again as described above.

Note: To abort calibration, you can press **meas** at any time.

Conductivity Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) The "TAG" menu allows the sensor TAG to be edited.
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 5) Perform the selected calibration as described on the following pages.
Follow the instructions on the display.

Note: Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

Auto Calibration

(Automatic calibration with specification of the calibration solution used)

NOTICE!

- Make sure that the values of the calibration solutions used correspond exactly to those specified in this manual.
If not, the resulting cell constant will be incorrect.
 - When calibrating in a liquid, make sure that the sensor, the separate temperature probe (if present) and the calibration solution have the same temperature.
Only this ensures that the cell constant is determined correctly.
-

- 1) Select the calibration solution:
 - **NaCl 0.01 mol/l**
 - NaCl 0.1 mol/l
 - NaCl sat.
 - KCl 0.01 mol/l
 - KCl 0.1 mol/l
 - KCl 1 mol/l
- 2) Press the **Start** softkey.
- 3) Immerse the sensor in the solution and press **Continue**.
- 4) Finally, the calibration data record will be displayed.
You can **Apply** or **Discard** these values.

Note: To abort calibration, you can press **meas** at any time.

Cond

“Entry of Solution” Calibration

(Calibration by entering the conductivity with display of the cell constant)

- 1) Press the **Start** softkey.
- 2) Immerse the sensor in the solution.
- 3) Enter the temperature-corrected conductivity value and press **enter**.
- 4) Finally, you can **Apply** these values or **Cancel** the calibration.

Cell Constant / Cell Factor Calibration

(Calibration by entering the cell constant (cell factor) with display of the conductivity)

- 1) Press the **Start** softkey.
- 2) Immerse the sensor in the solution.
- 3) Modify the value of the cell factor (cell constant) until the temperature-corrected conductivity value is reached. Then press **enter**.
- 4) Finally, you can **Apply** these values or **Cancel** the calibration.

Contacting conductivity sensor (conductive)	Cell constant
SE 202	0.100/cm \pm 2 %
SE 204	0.475/cm \pm 1,5 %
ZU 6985	1.19/cm \pm 1 %
SE 215 MS	1.00/cm \pm 2 %
Toroidal conductivity sensor (inductive)	Cell factor
SE 680 MS	6.4/cm

Installation Factor Calibration

- 1) Make sure that the sensor is in normal mounting position in the medium.
- 2) Press the **Start** softkey.
- 3) Modify the installation factor until the correct conductivity value is displayed (reference measurement). Then press **enter**.
- 4) Finally, you can **Apply** these values or **Cancel** the calibration.

Zero Calibration

- 1) Make sure that the sensor is outside the medium (in air).
- 2) Press the **Start** softkey.
- 3) Finally, you can **Apply** these values or **Cancel** the calibration.

Note: To abort calibration, you can press **meas** at any time.

Oxygen Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the desired "Calibration mode" and confirm by pressing **enter**.
- 4) Select "Membrane module replacement" if you wish to save a change of membrane or electrolyte in the connected sensor. The digital, optical oxygen sensor automatically recognizes when its sensor cap has been replaced.
- 5) The "TAG" menu allows the sensor TAG to be edited.
To do so, set "TAG" to **On** in the configuration menu (default setting: **Off**).
- 6) Perform the selected calibration as described on the following pages.
Follow the instructions on the display.

Note: Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

Calibration in Air

(Calibrating the slope in air)

- 1) Place sensor in air and wait for a stable measured value.
- 2) Press the **Start** softkey.
- 3) Adjust the correct value for "Relative humidity". Then press **Continue**.
Calibration will be performed.
- 4) Finally, you can **Apply** or **Discard** these values.

Note: To abort calibration, you can press **meas** at any time.

Oxy

Zero Calibration

(Zero calibration with oxygen-free medium, e.g., nitrogen 5.0)

- 1) Place sensor in oxygen-free medium and wait for a stable measured value.
- 2) Press the **Start** softkey. Calibration will be performed.
- 3) Finally, you can **Apply** these values or **Cancel** the calibration.

Data Entry Calibration

(Calibration by entering known sensor values)

- 1) Press the **Start** softkey.
- 2) Adjust the known sensor values for zero and slope.
- 3) Finally, you can **Apply** these values or **Cancel** the calibration.

Note: To abort calibration, you can press **meas** at any time.

pH

ORP

Oxy

Cond

Once you have completed all preparations, you can start with the actual measurement.

- 1) Connect the desired sensor to the meter. Some sensors require a special preparation. Please proceed according to the operating instructions for the sensor.
- 2) Switch the meter on using the **on/off** or **meas** key.
- 3) Depending on the measurement method and the sensor used, immerse the sensing part of the sensor in the medium to be measured.
- 4) Watch the display and wait for the reading to stabilize.

Note: Measurement can also be controlled via the Paraly SW 112 software.

Switching the Measured Value Display

During measurement, you can toggle between display of primary / secondary measured values and clock by pressing **meas**.

Adjusting the Temperature

When you connect a sensor without temperature detector, you can manually adjust the temperature for measurement or calibration:

- 1) Press **meas** to access measuring mode.
The adjusted temperature will be displayed.
- 2) Set the desired temperature value using the ▼ or ▲ arrow.
Holding the key depressed changes the temperature value at high speed.

pH

ORP

Oxy

Cond

The Data Logger

The meter provides a data logger. **Prior to use**, it must be configured and then activated. You can choose from the following logger types:


- Shot (manual logging by pressing the **Save value** softkey)
- Interval (time-controlled logging at a fixed interval)
- Difference (signal-controlled logging of measured variable and temperature)
- Intv+Diff (combined time- and signal-controlled logging)
- Limit value (combined time- and threshold-controlled logging)

The data logger records up to 10 000 entries, which can be assigned to different points of measurement (TAGs) and notes. The following data will be recorded: meas. point, note, sensor ID, serial number of sensor (Memosens), primary value, temperature, time stamp, device status.

The Paraly SW 112 software allows convenient management of the data logger.

It is always the currently selected process variable which is recorded.

Display: Icons related to the data logger

pH

ORP

Oxy

Cond


Operating Modes of the Data Logger (Logger Type)

Shot

In this mode, a measured value is recorded when the **Save value** softkey is pressed. In the measuring mode (**meas**), it is always possible to hold a value and then save it.


Interval (time-controlled)

In the "Interval" mode, the data are cyclically recorded.

Difference

When the delta range (process variable and/or temperature) related to the last entry is exceeded, a new entry is created and the delta range is displaced upwards or downwards by the delta value. The first entry is automatically created when the data logger is started.

pH


ORP

Oxy


Cond

Interval and difference (combined)

When the delta range related to the last DIFF entry is exceeded, a new entry is created (example: entry **A**) and the delta range is displaced upwards or downwards by the delta value. As long as the measured value remains within the delta range, logging is performed at the preset interval. The first DIFF entry is automatically created when the data logger is started.

**Limit value (combined)**

When one of the two limit values (Min/Max) is exceeded, the data are logged as defined by the "event interval". Additionally, the last ten measured values before an event are recorded (pre-trigger). As long as the measured value remains within the limits, logging is performed at the preset "basic interval".

pH

ORP

Oxy

Cond

Configuring the Data Logger

Prerequisite: Data logger is stopped.

The “Data logger” menu shows the number of occupied entries as well as the number of free entries. Configuration can also be done in the “Configuration” menu under “Data logger”.

- 1) Press **Menu** softkey.
- 2) Select “Data logger” and confirm by pressing **enter**.
- 3) Select “Configure data logger” and confirm by pressing **enter**.
- 4) Configure data logger as required (see table).
- 5) When you have completed the configuration, you can start the data logger!

Increasing the Battery Life

To increase the battery life for logger operation, the time for the display lighting selected in the configuration should be as short as possible.

Note: When the selected time has expired, display and backlighting switch off automatically. They can be switched on again by pressing any key.

pH

ORP

Oxy

Cond

Configuring the data logger (default in bold print)

Meas.point	Without		
Note	Without		
Recording	Non-circular		
	Circular		
Logger type	Shot		
	Interval	Interval	00:00:01...12:59:59 00:02:00
	Difference	1st difference	On Off
		Delta pH	pH 0.0...16.0 pH 1.0
		Delta mV	0 ... 2000 mV 1 mV
		Delta cond	0 ... 2000 mS/cm 1.0 µS/cm
		Delta conc	0 ... 9.99 % by wt 1%
		Delta MΩcm	0 ... 9.999 MΩcm 1.0 MΩcm
		Delta salinity	0.00 ... 45.0 g/kg 1.0 g/kg
		Delta TDS	0.00 ... 2000.0 mg/l 1 mg/l
		Delta saturation	0 ... 200% Air 1% Air
		Delta conc	0 ... 20.0 mg/l 1 mg/l
	Delta mbar	0 ... 1000 mbar 1 mbar	
		2nd difference	On Off
		Delta °C	0...99.9 °C 1.0 °C
		Delta °F	0...450 °F 1.0 °F
	Intv+Diff	Interval	see logger type: interval
Difference		see logger type: difference	
Limit value	Interval	Basis 00.00.01...12:59:59 00:01:00	
		Event 00.00.01 ...12:59:59	
	Limit values	Min/Max corresponding to permissible range (see Specifications)	

pH

ORP

Oxy

Cond

Starting/Stopping the Data Logger

With the data logger activated, automatic switch-off is disabled.

Every time the meter has been switched off, the data logger must be restarted.

Depending on the assignment of the right softkey (see "Configuration", page 14), you can start/stop the data logger as follows:

Right softkey	
Start/Stop logger	1) Press right softkey Start logger / Stop logger .
Hold value	2) Press Menu softkey. 3) Select "Data logger" using the arrows and confirm by pressing enter . 4) Press Start or Stop softkey, resp.

Viewing the Logger Data

In the "Data logger" menu you can view the recorded entries either individually or as curve characteristic (see examples).

You can also use the Paraly SW 112 software for managing the data logger.

- 1) Press **Menu** softkey.
- 2) Select "Data logger" using the arrows and confirm by pressing **enter**.
- 3) Select "View logger data" using the arrows and confirm by pressing **enter**.
- 4) Select filter ("Meas.point" or "Time + Meas.point" or "All values").
- 5) Select the parameter corresponding to the sensor.
- 6) Press **Menu** softkey.
- 7) Select the desired entries using the arrow (see example 1).
- 8) For display as curve characteristic, press **Graphic** softkey.

You can use the arrows to navigate between entries (see example 2).

Deleting the Logger Data

To delete the recorded entries, proceed as follows:

- 1) Press **Menu** softkey.
- 2) Select "Data logger" using the arrows and confirm by pressing **enter**.
- 3) Select "Delete logger data" using the arrows and confirm by pressing **enter**.
- 4) Select deletion mode: "Complete", "Data", "Meas.point" or "Filter"
(you can filter for measuring point, parameter or time).
- 5) Press **Delete** softkey. The data are deleted according to the configuration.
- 6) Press **Back** softkey to return to menu selection.

Example 1: Viewing the logger data

Data Logger

Primary measured value → **pH 7.000** **25.0 °C** ← Measured temperature

☺ **15.08.2012 10:45:36** ← Date and time

Meas.point ABC

Note

Sensor Analog

Scroll using arrow keys → **1/3**

Back **Graphic**

Example 2: Curve characteristic

Data Logger

7.50 25.8

pH °C

6.50 24.3

Entry 2 of 3 → **2/3**

Select entry using arrows. →

Back **Back** **Back**

Measured values Date Time

pH 7.21 25.0 °C 24/08/2012 17:12:50


pH

ORP


Oxy

Cond

Error and status messages appear as plain text on the display. By pressing **enter** and **Help**, more detailed help texts can be displayed. Information on the sensor condition is indicated by the “Sensoface” icon (friendly, neutral, sad) possibly accompanied by an info text.

Example of an error message: Press **enter** and **Help** to access the help text.

Help text for error 21

Sensoface (the “smiley” icon) provides information on the sensor condition (maintenance request). Measurement can still be performed. After a calibration, the corresponding Sensoface icon (friendly, neutral, sad) is shown together with the calibration data. Otherwise, Sensoface is only visible in measuring mode.

pH

ORP

Oxy

Cond

“Sensoface” Messages

The “Sensoface” icon provides information on the sensor condition:

Sensoface Meaning

Sensor is okay

Calibrate the sensor soon

Calibrate or replace the sensor

Info and Help Texts

When an error or status message appears on the screen, proceed as follows to view the corresponding info or help text:

- 1) Press **enter**.
- 2) Press the **Help** softkey.
- 3) The help text will be displayed. In most cases, you can remedy the cause of the error by yourself. Please refer to the following table for possible remedies.

Info	Message
Info 01	Cal timer expired
Info 02	Sensor wear
Info 03	Bad glass impedance
Info 05	Zero/Slope
Info 06	Response time too long
Info 07	Operating point (ISFET)
Info 08	Leakage current (ISFET)
Info 09	ORP offset
Info 10	Polarization

pH

ORP

Oxy

Cond

Error Messages

Error	Message	Remedy
 blinks	Replace the batteries	Replace the batteries.
ERR 1	Primary variable range	Check whether the measurement conditions correspond to the adjusted measuring range.
ERR 2	ORP range	
ERR 3	Temperature range	
ERR 4	Zero point	Thoroughly rinse the sensor and recalibrate.
ERR 5	Slope	If this does not help, replace the sensor.
ERR 6	Cell constant too high/ too low	Enter nominal cell constant or calibrate the sensor using a known solution.
ERR 7	Air pressure range	Check if the opening for the pressure sensor located on the back of the device is blocked.
ERR 8	Identical buffers!	Use a buffer solution with a different nominal value before starting the next calibration step.
ERR 10	Buffers interchanged!	Repeat calibration.
ERR 11	Unstable value (Drift too high)	Leave the sensor in the liquid until the measured value is stable. If this does not help, replace the sensor.
ERR 14	Time and date invalid	Set the date and time.
ERR 18	System error	Restart, reset to factory settings, configure and calibrate. If the error occurs again, contact the Service.
ERR 19	Factory settings error	Data error, measurement with analog sensors no longer possible. Contact the Service team.
ERR 21	No sensor connected	Connect operational Memosens sensor.
ERR 25	Buffer distance	Re-enter the buffer table (Paraly SW 112).
ERR 30	Data logger full	Clear the logger completely or partially.
ERR 31	MemoLog full	Clear the MemoLog completely or partially.

pH

ORP

Oxy

Cond

Option 001 SOP

The option meets specific requirements of the pharmaceutical and biotechnological industries.

SOP Cal

Here, you specify which buffers are to be used in which sequence.

You can combine buffer solutions from different buffer sets. Please note that the minimum distance allowed between two buffer solutions is $\Delta 2$ pH.

User Management (Access Control)

You can create up to 4 users with different access rights for configuration or calibration.

Sensor Verification

To make sure that only selected sensors can be operated on the meter, you can evaluate the sensor type and/or the "TAG" and "Group" data stored in the sensor.

The sensor will only be accepted if the data stored in the sensor corresponds to the data stored in the meter.

Temperature Adjustment

For Memosens sensors, you can perform a 1-point calibration of the internal temperature detector.

Option 002 Temp.Cal

(included in Option 001 SOP)

Temperature Adjustment

For Memosens sensors, you can perform a 1-point calibration of the internal temperature detector.

Enabling the Option 001 SOP

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Select Option "001 SOP" and enter your activation code.

Adapt SOP Cal

Here, you specify which buffers are to be used in which sequence.

You can combine buffer solutions from different buffer sets. Please note that the minimum distance allowed between two buffer solutions is $\Delta 2$ pH.

The "Configuration - Calibration" menu is extended as follows:


+ Calibration	
Cal mode	Calimatic Manual Data entry SOP cal
Adapt SOP cal	
Cal points	1-point 2-point 3-point
Buffer 1	
Buffer set	Mettler-Toledo 2.00/4.01/7.00/9.21 Knick CaliMat 2.00/4.00/7.00/9.00/12.00 Ciba 2.06/4.00/7.00/10.00 NIST technical 1.68/4.00/7.00/10.01/12.46 NIST standard 1.679/4.006/6.865/9.180 Hach 4.01/7.00/10.01/12.00 WTW 2.00/4.01/7.00/10.00 Hamilton 2.00/4.01/7.00/10.01/12.00 Reagecon 2.00/4.00/7.00/9.00/12.00 DIN 19267 1.09/4.65/6.79/9.23/12.75 Metrohm 4.00/7.00/9.00 User buffer 1
Buffer	Select a buffer from the selected set
Buffer 2	Select buffer set 2 and buffer (see buffer 1)
Buffer 3	Select buffer set 3 and buffer (see buffer 1)
Verification	Off On
Delta pH	pH 0.05 (Enter maximum permitted deviation from verification buffer; exceeding this value generates an error message)
Verification buffer	Select buffer set and buffer (see buffer 1)

Selecting SOP Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the "SOP cal" calibration mode and confirm by pressing **enter**.

Performing an SOP Calibration

The sequence of buffers to be used is displayed as specified in the configuration. After each calibration step, the identified buffer is marked off in the display. The next operation procedure is displayed. Perform the calibration following the instructions given in the display.

Note: Calibration is not possible when the device is connected via USB with the Paraly SW 112 software.

pH

ORP

Oxy

Cond

User Management (Access Control)

You can create up to 4 users with different access rights for configuration or calibration.

Enabling the User Management

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Select "User management"/"Management – Enable".
- 4) Select
User = ADMIN
PIN code = 1989 (factory setting)
- 5) Press **enter**.
- 6) To create more users / assign more PIN codes: Press **Continue** softkey.

Creating a User / Changing the PIN Code

1. Select a user (e.g., "User 1", default: ADMIN, PIN code 1989):

Configuration	
Management – Disable	
+ User 1	
+ User 2	
+ User 3	
+ User 4	
Return	Continue

Configuration	
Management – Disable	
- User 1	
ADMIN	
PIN code	1989
cal level	Access
conf level	Access
Return	Continue

You can allow or block access to configuration or calibration for each user.

2. Select ADMIN to open the editor for entering the user name:

Configuration	
Management – Disable	
- User 1	
ADMIN	
PIN code	1989
cal level	Access
conf level	Access
Return	Continue

Configuration											
End											
ADMIN											
1	2	3	4	5	6	7	8	9	0	?	`
Q	W	E	R	T	Z	U	I	O	P		*
A	S	D	F	G	H	J	K	L			'
>	Y	X	C	V	B	N	M	:	:	_	-
abc	Finish										

NOTICE!

If you lose the PIN code for the ADMIN user, the system access will be blocked. The manufacturer can generate a rescue PIN code.

pH

ORP

Oxy

Cond

How to Enter the Rescue PIN Code

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" using the cursor keys.
- 3) Press the ▼ and ▲ keys simultaneously.
- 4) Set User to "ADMIN".
- 5) Select "PIN code" and enter the rescue PIN code.
Confirm by pressing **enter**.
- 6) Press the **Continue** softkey.

Sensor Verification

To make sure that only selected sensors can be operated on the meter, you can evaluate the following data stored in the sensor:

- Model (sensor model)
- TAG (e.g., point of measurement)
- Group (e.g., facility)

With Option 001 enabled, the "Configuration" menu is extended as follows:

- Sensor verification	
Check model	Off Info Reject
Check TAG	Off Info Reject
Check group	Off Info Reject

You can select the following options:

- Off** No verification.
- Info** When a wrong sensor is connected, an error message will be displayed.
Nevertheless, you can continue working with the sensor.
- Reject** Here you specify values with which the sensor will be rejected.

pH

ORP

Oxy

Cond

Enabling the Option 002 Temp.Cal (included in Option 001 SOP)

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Configuration" and confirm by pressing **enter**.
- 3) Select Option "002 Temp.cal" and enter your activation code.

Selecting Temperature Calibration

- 1) In measuring mode, press the **Menu** softkey.
- 2) Select "Calibration" and confirm by pressing **enter**.
- 3) Select the "Temperature" calibration mode and confirm by pressing **enter**.

For Memosens sensors, you can perform a 1-point calibration of the internal temperature detector. To do so, enter the reference temperature and confirm the temperature adjustment by pressing the **Apply** softkey:

The screenshot shows a 'Calibration' screen with a red header. The current temperature is 24.4 °C. Below this, it prompts the user to 'Enter the reference temperature.' The 'Temp. offset' is set to -0.3 K. The 'Reference temperature' is set to 24.1 °C, which is highlighted with a red box. At the bottom, there are two buttons: 'Cancel' and 'Apply'.

Calibration	
24.4 °C	
Enter the reference temperature.	
Temp. offset	-0.3 K
Reference temperature	24.1 °C
Cancel	Apply

pH

ORP

Oxy

Cond

Connections	2x socket, 4 mm dia., for separate temp. detector 1 x M8 socket, 4 pins, for Memosens lab cable 1x micro USB-B for data transmission to PC 1 x socket depending on device version: Portavo 907 MULTI PH: pH socket acc. to DIN 19 262 Portavo 907 MULTI COND: Multi-contact for 2-/4-el. sensors Portavo 907 MULTI OXY: M12, 8 pins, for Memosens sensors or SE 340 sensor (optical oxygen)	
Air pressure measurement	700 ... 1100 hPa	
User interface	Straightforward menu navigation with graphic icons and detailed operating instructions in plain text	
Languages	German, English, French, Spanish, Italian, Portuguese, Russian	
Sensoface	Status indication (friendly, neutral, sad)	
Status indicators	For battery power level, logger	
Graphic display	QVGA TFT display with white backlighting	
Keypad	[on/off], [meas], [enter], [◀], [▶], [▲], [▼], 2 context-sensitive softkeys	
Data logger	10 000 memory locations	
Recording	Manual, interval- or event-controlled, with management of tag numbers and notes	
Cal data logger MemoLog (Memosens only)	Up to 100 Memosens calibration records can be saved	
	Recording	Directly retrievable via MemoSuite or Paraly SW 112 (USB)
	Viewable on the display	Manufacturer, sensor type, serial no., zero, slope, calibration date
Temperature input	2 x 4 mm dia. for integrated or separate temperature detector	
Measuring ranges	NTC30 temperature detector	-20 ... +120 °C / -4 ... +248 °F
	Pt1000 temperature detector	-40 ... +250 °C / -40 ... +482 °F
Measuring cycle	Approx. 1 s	
Measurement error ^{1,2,3)}	< 0.2 K (Tamb = +23 °C / +73.4 °F); TC < 25 ppm/K	

1) according to EN 60746-1, at nominal operating conditions

2) ± 1 count

3) plus sensor error

pH

ORP

Oxy

Cond

Communication	USB 2.0
Profile	HID, driverless installation
Usage	Data exchange and configuration via Paraly SW 112 software

Diagnostics functions

Sensor data (Memosens only)	Manufacturer, sensor type, serial number, wear, operating time, remaining lifetime, max. temperature, adaptive calibration timer, calibration and adjustment data, SIP, CIP and autoclaving counter
Calibration data	Calibration date; pH/Oxy: zero, slope; Cond: Cell constant
Device self-test	Automatic memory test (FLASH, EEPROM, RAM)
Device data	Device type, software version, hardware version

Data retention	Parameters, calibration data > 10 years
-----------------------	---

EMC	EN 61326-1 (General Requirements)
Emitted interference	Class B (residential environment)
Immunity to interference	Industrial environment EN 61326-2-3 (Particular Requirements for Transmitters)

RoHS conformity	According to directive 2011/65/EU
------------------------	-----------------------------------

Power supply	4 x AA alkaline batteries or 1x Li-ion battery, USB chargeable
---------------------	---

Nominal operating conditions

Ambient temperature	-10 ... +55 °C / +14 ... +131 °F
Transport/ Storage temperature	-25 ... +70 °C / -13 ... +158 °F
Relative humidity	0 ... 95 %, short-term condensing allowed

Housing

Material	PA12 GF30 (silver gray RAL 7001) + TPE (black)
Protection	IP66/67 with pressure compensation
Dimensions	Approx. (132 x 156 x 30) mm
Weight	Approx. 500 g

pH

ORP

Analog pH/mV input	pH socket, DIN 19 262 (13/4 mm)	
pH range	-2 ... 16	
Decimal places *)	2 or 3	
	Input resistance	1 x 10 ¹² Ω (0 ... 35 °C)
	Input current	1 x 10 ⁻¹² A (at RT, doubles every 10 K)
Measuring cycle	Approx. 1 s	
Measurement error ^{1,2,3)}	< 0.01 pH, TC < 0.001 pH/K	
mV range	-1300 ... +1300 mV	
Measuring cycle	Approx. 1 s	
Measurement error ^{1,2,3)}	< 0.1 % meas. val. + 0.3 mV, TC < 0.03 mV/K	
Memosens pH input (also ISFET)	M8 socket, 4 pins, for Memosens lab cable or M12 socket for Memosens sensors (Portavo 907 MULTI OXY only)	
Display ranges ⁴⁾	pH	-2.00 ... +16.00
	mV	-1999 ... +1999 mV
	Temperature	-50 ... +250 °C / -58 ... +482 °F
Sensor standardization *	pH calibration	
Operating modes *	Calimatic	Calibration with automatic buffer recognition
	Manual	Manual calibration with entry of individual buffer values
	Data entry	Data entry of zero and slope
Calimatic buffer sets *	-01- Mettler-Toledo	2.00/4.01/7.00/9.21
	-02- Knick CaliMat	2.00/4.00/7.00/9.00/12.00
	-03- Ciba (94)	2.06/4.00/7.00/10.00
	-04- NIST technical	1.68/4.00/7.00/10.01/12.46
	-05- NIST standard	1.679/4.006/6.865/9.180
	-06- HACH	4.01/7.00/10.01/12.00
	-07- WTW techn. buffers	2.00/4.01/7.00/10.00
	-08- Hamilton	2.00/4.01/7.00/10.01/12.00
	-09- Reagecon	2.00/4.00/7.00/9.00/12.00
	-10- DIN 19267	1.09/4.65/6.79/9.23/12.75
	-11- Metrohm	4.00/7.00/9.00
	-U1- (User)	loadable via Paraly SW 112
Permissible calibration range	Zero point	6 ... 8 pH
	With ISFET:	-750 ... +750 mV
	Operating point (asymmetry)	
	Slope	Approx. 74 ... 104 %
	(possibly restricting notes from Sensoface)	
Calibration timer *	Interval 1 ... 99 days, can be switched off	
Sensoface	Provides information on the sensor condition	
Evaluation of	zero/slope, response, calibration interval	

* user-defined

1) according to EN 60746-1, at nominal operating conditions

2) ± 1 count

3) plus sensor error

4) ranges depending on Memosens sensor

Memosens input	M8 socket, 4 pins, for Memosens lab cable or	
ORP	M12 socket for Memosens sensors (Portavo 907 MULTI OXY only)	
Display ranges ⁴⁾	mV	-1999 ... +1999 mV
	Temperature	-50 ... +250 °C / -58 ... +482 °F
Sensor standardization *	ORP calibration (zero adjustment)	
Permissible calibration range	Δ mV (offset)	-700 ... +700 mV

* user-defined

4) Ranges depending on Memosens sensor

Cond

Conductivity input, analog	Multi-contact for 2-/4-electrode sensors with integrated temp detector	
Measuring ranges	SE 202 sensor	0.01 ... 200 µS/cm
	SE 204 sensor	0.05 ... 500 mS/cm
	2-electrode sensors	0.1 µS * c ... 200 mS * c ⁴⁾
	4-electrode sensors	0.1 µS * c ... 1000 mS * c ⁴⁾
Permissible cell constant	0.005 ... 200.0 cm ⁻¹ (adjustable)	
Measurement error ^{1,2,3)}	< 0.5 % meas.val. + 0.4 µS * c ⁴⁾	
Conductivity input, Memosens	M8 socket, 4 pins, for Memosens lab cable or M12 socket for Memosens sensors (Portavo 907 MULTI OXY only)	
Measuring range	SE 615/1-MS sensor	10 µS/cm ... 20 mS/cm
Conductivity inputs		
Measuring cycle	Approx. 1 s	
Temperature compensation	Linear 0 ... 20 %/K, reference temperature adjustable nLF: 0 ... +120 °C / +32 ... +248 °F NaCl (ultrapure water with traces) HCl (ultrapure water with traces) NH ₃ (ultrapure water with traces) NaOH (ultrapure water with traces)	
Display resolution (autoranging)	Conductivity	0.001 µS/cm (c < 0.05 cm ⁻¹) 0.01 µS/cm (c = 0.05 ... 0.2 cm ⁻¹) 0.1 µS/cm (c > 0.2 cm ⁻¹)
	Resistivity	00.00 ... 99.99 MΩ cm
	Salinity	0.0 ... 45.0 g/kg (0 ... +30 °C / +32 ... +86 °F)
	TDS	0 ... 1999 mg/l (+10 ... +40 °C / +50 ... +104 °F)
	Concentration	0.00 ... 100 % by wt
Concentration determination	NaCl	0 – 26 wt% (0 °C / +32 °F) ... 0 – 28 wt% (+100 °C / +212 °F)
	HCl	0 – 18 wt% (-20 °C / -4 °F) ... 0 – 18 wt% (+50 °C / +122 °F)
	NaOH	0 – 13 wt% (0 °C / +32 °F) ... 0 – 24 wt% (+100 °C / +212 °F)
	H ₂ SO ₄	0 – 26 wt% (-17 °C / -1.4 °F) ... 0 – 37 wt% (+110 °C / +230 °F)
	HNO ₃	0 – 30 wt% (-20 °C / -4 °F) ... 0 – 30 wt% (+50 °C / +122 °F)
	H ₂ SO ₄	94 – 99 wt% (-17 °C / -1.4 °F) ... 89 – 99 wt% (+115 °C / +239 °F)
	HCl	22 – 39 wt% (-20 °C / -4 °F) ... 22 – 39 wt% (+50 °C / +122 °F)
	HNO ₃	35 – 96 wt% (-20 °C / -4 °F) ... 35 – 96 wt% (+50 °C / +122 °F)
	H ₂ SO ₄	28 – 88 wt% (-17 °C / -1.4 °F) ... 39 – 88 wt% (+115 °C / +239 °F)
	NaOH	15 – 50 wt% (0 °C / +32 °F) ... 35 – 50 wt% (+100 °C / +212 °F)
Sensor standardization	Cell constant	Input of cell constant with simultaneous display of conductivity value and temperature
	Input of solution	Input of conductivity of the calibration solution with simultaneous display of cell constant and temperature
	Auto	Automatic determination of the cell constant with KCl solution or NaCl solution

1) acc. to EN 60746-1, at nominal operating conditions

2) ± 1 count

3) plus sensor error

4) c = cell constant

Memosens input, amperometric oxygen	M8 socket, 4 pins, for Memosens lab cable or M12 socket for Memosens sensors (Portavo 907 MULTI OXY only)
Display ranges ⁴⁾	Saturation 0.000 ... 200.0 % Concentration 000 µg/l ... 20.00 mg/l Partial pressure 0.0 ... 1000 mbar
Temperature meas. range ⁴⁾	-20 ... +150 °C / -4 ... 302 °F
Sensor standardization	Automatic calibration in air (100 % RH) Zero calibration
Storage	in quiver with moistening sponge
Input optical oxygen	M12 socket for SE 340 sensor or Memosens sensors
OXY meas. ranges at +20 °C / +68 °F	Saturation 0.000 ... 200.0 % Concentration 000 µg/l ... 20.00 mg/l Partial pressure 0.0 ... 1000 mbar
Response time	t ₉₀ < 30 s t ₉₉ < 60 s
Measurement error ^{1,2,3)}	Zero signal < 0.1 % full saturation
Temperature meas. range ⁴⁾	0 ... +50 °C / +32 ... +122 °F
Measurement error ^{1,2,3)}	Temperature ± 0.2 K
Sensor standardization	Automatic calibration in air Zero calibration
Max. overpressure	2.5 bar
Immersion depth	Min. 60 mm Max. 25 m
Storage	in quiver with moistening sponge

1) according to EN 60746-1, at nominal operating conditions

2) ± 1 count

3) plus sensor error

4) ranges depending on Memosens sensor

A

AA batteries 6
 Access control (Option 001 SOP) 45
 Adapt SOP cal 43
 ADMIN PIN code lost 45
 Arrow keys 9
 Auto, conductivity calibration 27
 Automatic calibration (conductivity) 27
 Automatic calibration (pH) 23
 Auto, pH calibration 23

B

Basic interval (limit value logger type) 34
 Battery capacity 6
 Battery compartment 6
 Battery icon 6
 Battery life, increasing 35
 Battery replacement 6
 Buffer set (pH configuration) 14

C

Calibration (Cond), auto 27
 Calibration (Cond), cell constant 28
 Calibration (Cond), entry of solution 28
 Calibration, conductivity 27
 Calibration in air (Oxy) 29
 Calibration, ISFET 26
 Calibration, ORP 25
 Calibration (Oxy), data input 30
 Calibration, oxygen 29
 Calibration (Oxy), zero point 30
 Calibration (pH), Calimatic 23
 Calibration (pH), data input 24
 Calibration (pH), manual 24
 Calibration, pH/ORP combo sensor 24
 Calibration points (pH) 23
 Calibration record 10
 Calibration, SOP Cal (Option) 42
 Calibration, Temp.cal (Option) 47
 Calibration, zero point (operating point), ISFET 26
 Calimatic calibration 23
 Cal mode 26
 Cell constant, conductivity calibration 28
 Charge level of batteries 6
 CIP (sensor information) 10
 Clearing the logger 37
 Combo sensor, pH/ORP, calibration 24
 Conductivity calibration 27
 Conductivity configuration 18
 Configuration (conductivity) 18
 Configuration (ORP) 16
 Configuration (oxygen) 21
 Configuration (pH) 14
 Configuring the data logger 35
 Connecting a sensor 7
 Connecting cable for Memosens 7
 Connections 7
 Connection, USB (battery) 6
 Creating a user (Option 001 SOP) 45
 Curve characteristic (data logger) 38

D

Data input, Oxy calibration 30
 Data input, pH calibration 24
 Data logger configuration 35
 Data logger icons 32
 Data logger, operating modes 33
 Data logger, starting 37
 Data logger, stopping 37
 Data of the meter 48
 Deleting the logger data 37
 Delta range (data logger) 33
 Device info (Information menu) 13
 Device self-test 13
 Device test (Information menu) 13
 Difference (logger type) 33
 Display 9
 Display icons 8
 Displaying the logger data 37
 Disposal 5

E

Electrolyte, Oxy calibration 29
 Enabling an option 43
 Entry of solution, conductivity calibration 28
 ERROR (error codes) 41
 Error messages, overview 41
 Event interval (limit value logger type) 34

H

Help texts 40

I

Icons 8
 Icons for data logger 32
 In air, Oxy calibration 29
 Increasing the battery life 35
 Information (menu) 10
 Info texts 40
 Inserting the batteries 6
 Installation factor, calibration 28
 Interval and difference (logger type) 34
 Interval (logger type) 33
 ISFET calibration 26

K

Keypad 9

L

Limit value (logger type) 34
 Lithium-ion battery (start-up) 6
 Logger type "Difference" 33
 Logger type "Interval" 33
 Logger type "Interval and difference" 34
 Logger type "Limit value" 34
 Logger type "Shot" 33
 Loss of PIN code 45

M

Manual calibration (pH) 24
meas, switch-on 8
Measuring 31
Membrane module replacement 29
MemoLog (Memosens only) 12
Memosens connecting cable 7
Memosens sensors 7
Messages (Information menu) 12
Micro USB port 7

N

Note (data logger) 36

O

on/off, switch-on 8
Operating modes of data logger 33
Operating time of cap (optical oxygen) 10
Operating time of sensor (sensor information) 10
Option 001 SOP 42
Option 001 SOP, enabling 43
Option 002 Temp.cal 42
Option 002 Temp.cal, enabling 47
ORP calibration 25
ORP configuration 16
Overview of error messages 41
Overview of status messages 39
Oxygen calibration (Oxy) 29
Oxygen configuration 21

P

Passcode lost 45
pH calibration 23
pH configuration 14
pH/ORP sensor, calibration 24
PIN code, changing (Option 001 SOP) 45
PIN code lost 45
Ports 7
Power-on 8
Pressure correction (Oxy configuration) 21
Pre-trigger (limit value logger type) 34

R

Rechargeable battery, Li-ion 6
Redox calibration 25
Registered trademarks 5
Replacing the electrolyte (Oxy) 29
Replacing the membrane (Oxy) 29
Rescue PIN code entry 46

S

Saving a membrane module replacement (Oxy) 29
Sensoface messages 40
Sensor connection 7
Sensor information 10
Sensor monitor 12
Sensor network diagram 11
Sensor verification (Option 001 SOP) 46

Sensor without temperature detector 31
Serial number of cap (optical oxygen) 10
Serial number of sensor (sensor information) 10
Setup, conductivity 18
Setup, ORP 16
Setup, oxygen 21
Setup, pH 14
Shot (logger type) 33
SIP (sensor information) 10
Slope calibration, oxygen 29
Softkey 9
SOP calibration (Option 001) 44
Specifications 48
Status messages, overview 39
Stopping the data logger 37
Switching between displays 31
Switching on the meter 8
Symbols in display 8

T

Table of error messages 41
Table of info texts 40
TAG entry (Cond) 27
TAG entry (ISFET) 26
TAG entry (ORP) 25
TAG entry (Oxy) 29
TAG entry (pH) 23
Technical data 48
Temp.cal (Option) 47
Temperature, manual adjustment 31
Temperature probe, connection 7
Toggling the measured value display 31
Trademarks 5

U

USB port (battery) 6
USB port, micro 7
User management (Option 001 SOP) 45

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The screenshot shows the main menu of the 'Datenlogger' application. At the top, the title 'Datenlogger' is in a green bar. Below it, the main data is displayed: 'pH 7.000' and '25.0 °C'. A smiley face icon is to the left of the date and time '15.08.2012 10:45:36', which is circled in red. Below this, there are two columns of text: 'Messstelle ABC' and 'Notiz Sensor Analog'. At the bottom left, it says '1/3' with the text 'Mit Pfeiltasten blättern' pointing to it. At the bottom, there are two green buttons: 'Zurück' and 'Grafik'. Red arrows point from external text labels to various elements: 'Gemessene Hauptgröße' to 'pH 7.000', 'Gemessene Temperatur' to '25.0 °C', 'Datum und Uhrzeit' to the circled date and time, and 'Mit Pfeiltasten blättern' to '1/3'.

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The screenshot shows the 'Datenlogger' application in graph view. The title 'Datenlogger' is in a green bar. Below it is a line graph with two data series: 'pH' (red line) and '°C' (blue line). The y-axis has labels '7.50' and '6.50' on the left, and '25.8' and '24.3' on the right. A red oval highlights a specific data point on the graph, with an arrow pointing to it from the text 'Eintrag mit Pfeiltasten auswählen'. Below the graph, the selected entry's data is shown: 'pH 7.21 25.0 °C 24.08.2012 17:12:50'. At the bottom left, it says '2/3' with the text 'Eintrag 2 von 3' pointing to it. At the bottom, there are two green buttons: 'Zurück' and 'Grafik'. Red arrows point from external text labels to '2/3', 'pH 7.21 25.0 °C 24.08.2012 17:12:50', and the 'Zurück' button. The labels 'Messwerte', 'Datum', and 'Uhrzeit' are positioned below the data string with arrows pointing to the corresponding parts.

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