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# SE555 Series pH Sensors: Instructions for Use

**WARNING – Failure to observe this warning may result in serious injury.**

The safety alert symbol on the nameplate means:

**Read these Instructions for Use, observe the specifications, and follow the safety instructions.**

## 1 Safety Instructions

### 1.1 All Applications – All Sensor Variants

Hazards due to pressure, temperature, aggressive media or explosive atmospheres are possible, depending on the location of use. Installation, operation, and servicing of the sensor may therefore only be carried out by suitably trained personnel authorized by the operating company.

### 1.2 Hazardous Areas – All Sensor Variants

Observe all applicable local and national codes and standards for the installation of electrical equipment in hazardous locations. For orientation, please refer to IEC 60079-14, EU Directives 2014/34/EU and 1999/92/EC (ATEX), NFPA 70 (NEC), ANSI/ISA-RP12.06.01. The electrical and thermal parameters of the sensors must be adhered to.

### 1.3 Hazardous Areas – Sensors with Memosens Connector

Memosens Ex sensors are marked by an orange-red ring. Combined with a model CA/MS-\*\*\*X\*\* or model CA/MS-\*\*\*X\*\*-L measuring cable or a certified measuring cable which is identical in hardware and function, the sensor may be connected to a suitable measuring device, as described in the Certificates BVS 15 ATEX E141 X and IECEx BVS 15.0114X.

### 1.4 Hazardous Areas – Sensors with VP Connector

The sensors may only be operated with an approved device in floating, intrinsically safe circuits.

## 2 Intended Use

The sensor is used for continuous measurement of pH and, optionally, ORP (AMSN) in liquid media. The SE555 is a low-maintenance sensor with pressurized electrolyte and integrated temperature detector for automatic temperature compensation. The sensor is sterilizable by autoclaving and is CIP-/SIP-capable. The sensor is designed for applications in industrial processes:

- hygienic processes
- biotechnology, food, pharmaceuticals
- high temperatures, high pH values, galvanic processes

## 3 Installation and Commissioning

- On unpacking, check the sensor for mechanical damage. Report any damage to your Knick service team.
- Take off the watering cap and use the included knife to remove the silicone sealing from the junction. Briefly rinse the sensor with pure water. After rinsing, the sensor should only be dabbed dry. Do not rub the pH-sensitive glass, since this can lead to electrostatic charging and sluggish response times.
- Check the space behind the pH-sensitive glass for the presence of any air bubbles and remove them by gently shaking the sensor up and down.
- Install the sensor in the fitting as specified in the User Manual for the respective fitting.
- Connect the sensor to the cable.

## 4 Operation

### 4.1 Calibrating the Sensor

2-point calibration is recommended for the SE555. First remove the watering cap. Then dip the sensor successively into two different buffer solutions with given pH values (e.g., CaliMat pH 7.00 and pH 4.00) and calibrate the sensor to these buffer values. ORP calibration is carried out as a 1-point calibration (delta value) with an ORP calibration solution. Please refer to the user manual of the pH transmitter for further details.

### 4.2 Sterilizing

For application in sterile processes, such as fermentation, sterilize the sensor before starting the operating cycle. Sterilization can be effected in situ with steam or superheated process medium.

### 4.3 Temperature Detector

The integrated temperature detector is intended for automatic compensation of the pH signal and not for high-precision and safe temperature indication or control of the process temperature.

## 5 Maintenance and Cleaning

The sensor tip and junction must be carefully rinsed with pure water after each operating cycle. Under no circumstances must process medium be allowed to dry on these parts.

When the sensor is not in operation, store it, with the sensor tip and junction, well submerged in a container with electrolyte (3 mol/l KCl). If a sensor is stored dry for a few days by mistake, let it soak in electrolyte for several hours before use.

Contamination can be removed as follows:

Protein:	Pepsin/HCl solution
Lime:	Diluted acids
Silver sulfide:	Thiourea/HCl solution
Grease:	Surfactant solution/base

## 6 Specifications

### Model Code

The markings on each sensor or on the packaging label include the following information:

SE555X/1-NMSN	<b>Model designation</b>
1	<b>Connector</b>
	MS: Memosens
	VP: VarioPin
	S8: coaxial
	<b>Pt electrode (ORP)</b>
	A: with
	N: without
	<b>Length</b>
	1: 120 mm
	2: 225 mm
	3: 325 mm
	4: 425 mm
	<b>Ex approval</b>
	X: yes
	N: no

### Further Data

Measuring range for pH / ORP	0 ... 14 / ±1500 mV
Temperature	0 ... 135 °C
Relative process pressure	-1 ... 6 bar
Pressure resistance	12 bar
Junction	1x ceramic
Electrolyte	Viscous gel with internal pressure
Reference system	Ag/AgCl with silver ion trap
Sensor material	Omega glass High impedance for high-temperature applications, minimal alkali error, CIP-/SIP-capable
ORP sensor material	Platinum
Body material	Glass
Mounting	PG 13.5
Tightening torque	1 ... 3 Nm
Temperature detector	MS connector: NTC 30 kΩ VP connector: Pt1000 S8 connector: without


## 7 Disposal

Observe the applicable local or national regulations for disposal.

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TA-SE555X-KNXX04 20220809

# Hazardous Locations: Electrical and Thermal Parameters

## For Sensors with VP and S8 Connectors

**Certificate Number:**

PTB 14 ATEX 2004

**Marking:** II 1/2 G Ex ia IIC T6...T3 Ga/Gb**Electrical Parameters:**

- From a safety point of view, all circuits, including the shield wire and ground, are to be regarded as galvanically interconnected.
- Metallic process connection parts must be included in the local equipotential bonding system.
- The effective internal inductance  $L_i$  and capacitance  $C_i$  are negligibly low.

**For Sensors with Temperature Detector (VP Connector):**

Temperature class	Maximum input voltage $U_i$	Maximum input current $I_i$	Maximum input power $P_i$
T6	12 V	30 mA	50 mW
T4	15 V	80 mA	110 mW
T3	18 V	170 mA	200 mW

**For Sensors without Temperature Detector (S8 Connector):**

Temperature class	Maximum input voltage $U_i$	Maximum input current $I_i$	Maximum input power $P_i$
T6, T4, T3	18 V	170 mA	200 mW

**Thermal Parameters:**

Exceeding the standardized atmospheric conditions within the manufacturer's specifications, e.g., with regard to the ambient temperature and pressure, does not endanger the equipment's durability.

**For Sensors with Temperature Detector (VP Connector):**

The maximum heating of relevance with regard to the temperature class is uniformly based on the available power of the externally connected circuits and is 427 K/W in air. The maximum permissible ambient temperature is 40 °C.

Temperature class	Maximum input power $P_i$	Maximum permissible process temperature
T6	50 mW	55 °C
T4	110 mW	80 °C
T3	200 mW 110 mW	100 °C 135 °C

**For Sensors without Temperature Detector (S8 Connector):**

The maximum heating of relevance with regard to the temperature class is negligibly low. The maximum permissible ambient temperature is 60 °C.

Temperature class	Maximum input power $P_i$	Maximum permissible process temperature
T6	200 mW	80 °C
T4	200 mW	130 °C
T3	200 mW	135 °C

**Special Conditions**

None

## For Sensors with Memosens Connector

**Certificate Number:**BVS 16 ATEX E 037 X  
IECEX BVS 16.0030X**Marking:** II 1G  
Ex ia IIC T3/T4/T6 Ga**Thermal Parameters:**

Temperature class	Ambient temperature range $T_a$	Maximum permissible process temperature
T6	-20 °C < $T_a$ < +70 °C	70 °C
T4	-20 °C < $T_a$ < +120 °C	120 °C
T3	-20 °C < $T_a$ < +135 °C	135 °C

**Special Conditions**

- The cable and sensor may only be used within the ambient temperature range specified for the temperature class.
- The measuring cable, including its connecting head, must be protected from electrostatic charging if it passes through areas of Zone 0 (category 1G).
- The Memosens sensors may not be operated in electrostatically critical processing conditions. Intense vapor or dust flows directly impacting on the connection system must be avoided.
- Metallic process connection parts must be mounted at the installation site so that they are electrostatically conductive (< 1 MΩ).