

CaliMat

Neue pH-Pufferlösungen



- Gebrauchsfertige pH-Pufferlösungen mit weitem Wertebereich pH 2,00 / 4,00 / 7,00 / 9,00 / 12,00
- Rückführbar auf SRM / NIST
- Kostengünstig
- Sicher durch Farbcodierung
- Einsatz in Industrie-Prozessen und im Labor
- Praktische und gebrauchssichere Handhabung durch integrierten Kalibrierbecher
- Ideale Nachfüllmenge 3000 ml für automatische Kalibriersysteme mit Unical 9000
- Lieferung mit Quality Control Certificate

Knick Elektronische Messgeräte GmbH & Co. KG

Beuckestraße 22 · 14163 Berlin
Telefon: +49 30 801 91 - 0 · Fax: +49 30 801 91 - 200
knick@knick.de · www.knick.de

CaliMat

New pH Buffer Solutions

pH value (20 °C)	Volume	Order no.
2.00 ± 0.02	250 ml	CS-P0200/250
4.00 ± 0.02	250 ml	CS-P0400/250
	1000 ml	CS-P0400/1000
	3000 ml	CS-P0400/3000
7.00 ± 0.02	250 ml	CS-P0700/250
	1000 ml	CS-P0700/1000
	3000 ml	CS-P0700/3000
9.00 ± 0.02	250 ml	CS-P0900/250
	1000 ml	CS-P0900/1000
	3000 ml	CS-P0900/3000
12.00 ± 0.05	250 ml	CS-P1200/250
Buffer sets		
Set 4.00	3 x 250 ml	CS-PSET4
Set 7.00	3 x 250 ml	CS-PSET7
Set 9.00	3 x 250 ml	CS-PSET9
Set 4.00, 7.00, 9.00	je 250 ml	CS-PSET479



- Ready to use pH buffer solutions with a wide range of values pH 2.00 / 4.00 / 7.00 / 9.00 / 12.00
- Traceable to SRM / NIST
- Reasonable
- Sure through colorcoding
- For use in industrial process applications as well as in laboratories
- Convenient and safe handling thanks to integrated calibration beaker
- Ideal filling quantity of 3000 ml for automatic calibration systems with Unical 9000
- Supplied with quality control certificate